Used for the Airborne Disinfection of Surfaces after cleaning.

NFT 72-281 (V 2014) compliant

Phileas® *range compatible*

2 years shelf life

*Pharma
quality*

*No
residues*

*Ready
to use*

Phileasafe®

*Ready to use disinfectant solution, with a broad spectrum.
Stabilized hydrogen peroxide
Peracetic acid +/-1500 ppm*

ADVANTAGES :

- Stabilised by a physical process
- • 0,2 µ filtered (pharmaceutical quality)
- • • Food contact material
- • • • Proven efficiency for airborne disinfection according to norm NF T 72-281

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Tested with Phileas® 50D or Phileas 75

HOW TO USE : Ready to use product – Diffuse the solution with 12 ml/m³ (effective dose), preferably with PHILEAS® machines – Contact time : 2 to 3 h – Vent the room before using it again (20 times with new air changing).

PRECAUTIONS :

Do not use if allergic to any of the components. Do not apply to mucous membranes or wounds. If ingested, consult a GP and show him/her the label. Keep away from freezing temperatures. R36 irritant for the eyes. S2 keep away from children. S24/25/26 avoid contact with eyes, rinse with water if necessary. Unused products should be disposed of with household waste. The empty container should be rinsed and can be processed and recycled. The equipment should be cleaned with water. Residues should not be poured down the drains, nor in waterways. Keep away from light sources, between 10 and 30°C.

PACKAGING : RefPHISAF 10 : 10 L can
Palettization : 400 L in 10 L cans

RefPHISAF 6 : 6 x 1 L bottles

The use of a Major pump is recommended to transfer the solution from the can to a smaller recipient.