

MyPhileas®

Security and traceability



Devea, biodecontamination expert and developer of the MicroDrop Technology, offers an application system **MyPhileas**, available on tablet or on computer, to monitor remotely and ensure the complete operations **traceability**. Compliant with 21 CFR Part 11 requirements, **MyPhileas** is key in **user security** and is easy to navigate.

Key features

- Secure connection with Phileas® devices via private wifi, with up to 5 devices
- Creation of programs, naming of zones
- Launching disinfection cycles, with or without delay, cycle monitoring
- Generation and download of diffusion reports
- Complete and secure traceability
- Access restriction according to your specifications
- Secure, localized and unalterable data
- Audit trail
- 21 CFR Part 11 compliance

Why DEVEA?

- Phileas® / O2SAFE® technology with efficiency validated by the **Marketing Authorisation FR-2019-0071** and in compliance with the **EN 17-272 standard**
- 4-log₁₀, 5-log₁₀, 6-log₁₀ bioburden reduction according to specifications
- Biocide 100% biodegradable, residue-free, non-corrosive
- Micro-Drop Technology: creation of micro-droplets by centrifugation, for optimal dispersion and evaporation. Disinfection of all surfaces, even inaccessible
- Designed, manufactured and maintained in France

Easy disinfection control and traceability



1- Remote control

Private wifi between Phileas® devices and MyPhileas app

- Zone creation and modification
- Disinfection programming
- Disinfection cycle monitoring

2- Traceability and 21 CFR Part 11 compliance

Reporting and data integrity

- Secure access by identification and password, or connection via *Active Directory*
- Access to features restricted by user category
- *Audit trail*: archiving of each use or modification
- Impossible to falsify a report, to modify the time, etc.
- Reports are automatically saved and cannot be deleted

3- Access to maintenance data

Data history and update

- Access to machine information: serial number, hours of use, etc.
- Adjustment of several points on the machine
- Access to the history of operations by the maintenance team



P30461 - TDS MyPhileas -- October 2022 - EN

